

Invacare® TDX® SI Rehab Seat Power Wheelchair



Price List and Order Form

Effective May 18, 2009

For ease of ordering, contact customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone _____ P.O. # _____ Date _____

Ship to Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

CHAIR TYPE

- TDXSI TDX SI with MK6i™ (K0848)..... 5,970.00 ¹
- TDXSISEAT TDXSI Seat Only..... 755.00

Group 3 – K0848 Basic Package

- TDX SI with MK6i and Adjustable Rehab Seat
- 300 lb Weight Capacity
- 1431-3 14" X 3" Tires w/Foam Filled InsertsStandard
- 2 Pole MotorsStandard
- Flip Back Fixed Height ArmStandard
- MK660 60 Amp ControllerStandard
- SPJPLUS SPJ™+ Single Drive JoystickStandard
- 70NHD 70° Swingaway FootrestsStandard
- 8 Amp Off Board Battery ChargerStandard
- 1515 Push Button Style Seat Positioning StrapStandard

TRANSPORT TIE-DOWN

- TRRO Transport Ready Option..... 300.00
- TRBKTS Wheelchair Transport Brackets..... 250.00

*TRANSPORT TIE-DOWN

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lb crash dummy, which corresponds to a person with a weight of 114 lb to 209 lb.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*Transport Tie-Down option is not retrofittable to existing models and is not field serviceable.

TRRO FEATURES FOR TDX SI

Features Required for TRRO	Optional Features for TRRO
<ul style="list-style-type: none"> • TDX SI • 70° Swingaway, 70° Tapered, or Manual Center Mount Front Riggings • T-arms required with 12"-15" seat depth • 18F,19F,28F,29F for 12"-16" widths only available with 16" 17" or 18" deep • Adjustable Rehab Seat with Adjustable Seat Tilt & Back Angle • Sectional or Straight Rear Uprights or No Push Handles • 14" x 3" Tires with Foam Filled Inserts • 90° or 95° Back Angle • 18" – 24" ASBA Back Height • Any Seat Width • Any Seat Depth • 300 lb • Low Seat-to-Floor Height • 22NF Batteries • Any Color Frame Finish 	<ul style="list-style-type: none"> • Adjustable Angle Flip Up Footplates, Flip-Up Composite Footplates or Aluminum Footplates-Deluxe • Wheel Locks • Calf Strap • Heel Loops w/Ankle Straps • Impact Guards • One Sided Fork Package • SPJ+, MPJ+ Joystick • All PSR+ / PSF+ Joysticks • Proportional Attendant Control • "T" Handle Flexible Joystick Ext. • Left Hand Mounted Joystick • Straight Handle Flexible Joystick Ext. • Easy Remote Programmer • Reinforced Back Upholstery

USER WEIGHT LIMIT

- U300 User Weight <300 lbs **Standard**

WHEEL OPTION (PAIR)

- 1431-3 14" x 3" w/Foam Filled Inserts..... **Standard**

FORK OPTION

- FKPKG One-Sided Fork Package 500.00

Footnotes:

1. TDX SI comes standard with: 300 lb user weight capacity, 5 MPH top speed, 2-pole motor package and MK660 Controller.

Invacare® TDX® SI Rehab Seat Power Wheelchair

Accessories and Options

FRAME FINISH

CONTEMPORARY

- 24P Wet Black No Charge
- 60P Silver Metallic No Charge
- 71P Silver Vein No Charge
- 105P Midnight Blue No Charge
- 115P Black Prism No Charge

TOO HOT

- 30P Sunny Yellow No Charge
- 61P Electric Red No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge

WAY COOL

- 62P Electric Blue No Charge
- 104P Electric Teal No Charge
- 121P Lolly Pop Blue No Charge
- 122P Cosmic Blue No Charge
- 125P Emerald Green No Charge
- 127P Grasshopper Green No Charge

SEAT

SEAT-TO-FLOOR HEIGHT

- LOW 16.5 at 0° **Standard**
- MED 18.5 at 0° No Charge
- TALL 20.5 at 0° No Charge³

SEAT SELECTION

- ADJASBA Adj. Rehab Seat w/Adj. Seat Tilt & Back Angle..... **Standard**
- ADJRECL Manual Recline for Adjustable ASBA (E1226) . 744.00
- CURV Curved No Charge
- CTB Contoura Seating (E2620/E2621) 625.00
- OMCU Omit Contoura Seat Cushion..... Less 20.00⁴
- VAN Van Style Seat 70.00⁴

ADJUSTABLE SEAT WIDTH RANGE

- 1216W 12" to 16" Seat Width Range No Charge⁵
- 1620W 16" to 20" Seat Width Range No Charge
- 2022W 20" to 22" Seat Width Range No Charge

SEAT WIDTH⁶

- W12 Seat Width 12" 110.00⁷
- W13 Seat Width 13" 110.00⁷
- W14 Seat Width 14" 110.00⁷
- W15 Seat Width 15" No Charge⁷
- W16 Seat Width 16" No Charge
- W17 Seat Width 17" No Charge
- W18 Seat Width 18" No Charge
- W19 Seat Width 19" No Charge
- W20 Seat Width 20" No Charge
- W21 Seat Width 21" 329.00
- W22 Seat Width 22" 329.00
- VW16 Seat Width 16" No Charge⁸
- VW18 Seat Width 18" No Charge⁸
- VW20 Seat Width 20" 329.00⁸
- VW22 Seat Width 22" 329.00⁸

ADJUSTABLE SEAT DEPTH RANGE

- 1218D 12" to 18" Seat Depth Range No Charge⁵
- 1619D 16" to 19" Seat Depth Range No Charge
- 1922D 19" to 22" Seat Depth Range No Charge

FOOTNOTES

³ The long extension manual center mount front rigging to provide up to 17.5" KTH is available using specials #C-071025041441. The seat angle is set at 5 degrees.

⁴ Not available with Transport Tie Down options (TRRO). Seat -to- Floor Height 21".

⁵ Must select 12" – 18" Depth Range (1218D) when ordering the 12" – 16" Seat Width Range (1216W).

⁶ Contoura available in even Seat Width from 16" – 22".

⁷ User weight limit is 165 lb.

⁸ Seat Widths and Seat Depths for Van seat Semi Recline with Solid base.

SEAT DEPTH SETTING

- D12 Seat Depth 12" 110.00⁷
- D13 Seat Depth 13" 110.00⁷
- D14 Seat Depth 14" 110.00⁷
- D15 Seat Depth 15" No Charge⁷
- D16 Seat Depth 16" No Charge
- D17 Seat Depth 17" No Charge
- D18 Seat Depth 18" No Charge
- D19 Seat Depth 19" No Charge
- D20 Seat Depth 20" No Charge
- D21 Seat Depth 21" 529.00
- D22 Seat Depth 22" 529.00
- VD16 Seat Depth 16" No Charge⁸
- VD17 Seat Depth 17" No Charge⁸
- VD18 Seat Depth 18" No Charge⁸
- VD19 Seat Depth 19" No Charge⁸
- VD20 Seat Depth 20" 529.00⁸

VAN UPHOLSTERY

- GV Grey Vinyl No Charge⁹

SEAT STYLE

- SEMIRB Semi Recline with Solid Base No Charge¹⁰

FIXED ARM TYPE (PAIR)

- 18F Full Length Fixed Height Conv..... No Charge
- 28F Desk Length Fixed Height Conv..... **Standard**

ADJUSTABLE ARM TYPE (PAIR)

- 19F Full Length Adj Hgt Conv. (E0973) 230.00
- 29F Desk Length Adj Hgt Conv. (E0973)..... 230.00
- TADM T-arms Desk Lgth Adj Hgt Arm (8"-12")(E0973)230.00
- TAFM T-arms Desk Lgth Adj Hgt Arm (8"-12")(E0973)230.00
- VAH Full Length Adjustable Flip-Up (E0973) 230.00¹¹
- VFH Full Length Fixed Height **Standard**¹¹
- VAHD Desk Length Adjustable Height (E0973) 230.00¹¹
- VFHD Desk Length Fixed Height..... No Charge¹¹

MIXED HEIGHT ADJUSTABLE ARM TYPE (EACH)

- 19FR Full Length Adj Hgt Right (E0973) 115.00
- 29FR Desk Length Adj Hgt Right (E0973) 115.00
- 19FL Full Length Adj Hgt Left (E0973)..... 115.00
- 29FL Desk Length Adj Hgt Left (E0973)..... 115.00

CANTILEVER ARM TYPE (PAIR)¹²

- CLVF Full Length Arm.....\$350.00
- CLVD Desk Length Arm.....\$350.00

CANTILVER MIXED ARMTYPE (EACH)

- CLVFR Full Length Arm Right.....\$175.00
- CLVDR Desk Length Arm Right.....\$175.00
- CLVFL Full Length Arm Left.....\$175.00
- CLVDL Desk Length Arm Left.....\$175.00

ARMPAD OPTIONS (CANTILEVER ARMS)

- WF Waterfall Armpads..... **Standard**
- RCL Reclining Armpads..... No Charge

SEAT BACK TYPE

- SECT Sectional Rear Uprights..... **Standard**
- SRUP Straight Rear Uprights No Charge
- 9720P No Push Handles (Straight Canes) 82.00
- L8661 Omit Rear Uprights..... Less 20.00

SEAT BACK ANGLE DEGREE

- 90BA 90° **Standard**
- 95BA 95° No Charge

FOOTNOTES

⁹ Only available with VAN seat selection and SEMIRB seat style.

¹⁰ Only available with VAN seat selection and GV Grey Vinyl.

¹¹ Arms only available with VAN seat selection and SEMIRB seat style.

¹² Cantilever Arms only available with 300 lb. rated systems. Not available with BPO, Otto Bock or 12"-15" seat width. Must use hard back. See chart on last page for arm height adjustment range.

Invacare® TDX® SI Rehab Seat Power Wheelchair

Accessories and Options

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Pos. Strap50.00 ¹³
- 1515 Push Button Style Seat Pos. Strap..... **Standard**

CHEST POSITIONING STRAPS

- 7311BK Airline Buckle Chest Pos. Strap (E0960)..... 100.00 ¹²
- 9515 Push Button Style Chest Pos. Strap (E0960) 100.00
- 7321BK Hook & Loop Chest Pos. Strap (E0960)..... 100.00

BACK HEIGHT ¹⁴

- BH12 Back Height 12"No Charge
- BH13 Back Height 13"No Charge
- BH14 Back Height 14"No Charge
- BH15 Back Height 15"No Charge
- BH16 Back Height 16"No Charge
- BH17 Back Height 17"No Charge
- BH18 Back Height 18"No Charge
- BH19 Back Height 19"No Charge
- BH20 Back Height 20"No Charge
- BH21 Back Height 21"No Charge
- BH22 Back Height 22"No Charge
- BH23 Back Height 23"No Charge
- BH24 Back Height 24"No Charge

CONTOURA BACK HEIGHT

- TSBH25 Back Height 25"No Charge
- TSBH28 Back Height 28"No Charge

BACK HEIGHT RECLINER

- BHR18 Back Height 18"No Charge
- BHR20 Back Height 20"No Charge
- BHR22 Back Height 22"No Charge
- BHR24 Back Height 24"No Charge
- BHR26 Back Height 26"No Charge

ELECTRONICS

MK6™ CONTROLLERS

- MK660ACC Controller w/ Actuator Control **Standard** ¹⁵
- GTRAC G-Trac Module (K0108) 995.00

STANDARD DRIVER CONTROL

- SPJPLUS SPJ+ Single Drive Joystick **Standard**

EXPANDABLE CONTROLLER SYSTEMS

- PWH Harness for Expandable Electronics (E2313) 350.00 ¹⁶
- MPJM6U MPJ+ Multiple Drive Joystick (E2377) 750.00
- PSM6 Personal Joystick, Inductive on Front 1,130.00 ¹⁷
- PSRM6 Personal Joystick, Inductive on Rear 1,130.00 ¹⁷
- PBOD Push Button On/Off w/PB Drive Select.....No Charge ¹⁸
- PBSS Push Button On/Off w/Speed Select.....No Charge ^{18,19}
 - OR On/Off Mounted Right.....No Charge ¹⁸
 - OL On/Off Mounted Left.....No Charge ¹⁸
- ASLRDYM6 MK6 Display Only
(No Drive Control) (E2377) 750.00 ²⁰

**FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE FOLLOWING:
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.**

FOOTNOTES

- ¹³ 1311BK & 7311BK not available on 12" – 15" wide.
- ¹⁴ Back Heights 23" & 24" not available with 12"-16" widths. With junior seat widths 12"-16", the back height is adjustable.
- ¹⁵ G-Trac provides Improved Tracking Capabilities and Enhanced Traction. Does not require Special Motors. Not available w/ SPJ+ joysticks.
- ¹⁶ Harness is required to interface joystick or MK6i display with expandable controllers. Option is located under Electronics Accessories in configurator.
- ¹⁷ TOD-SS: Toggle On/Off-Drive Select, Speed Potentiometer Standard.
- ¹⁸ PSRM6 Only.
- ¹⁹ Drive Select through mode switch, or mode switch and Joystick Command.
- ²⁰ Display Only, No Driver Control. Allows choice from order form.

ALTERNATIVE PROPORTIONAL CONTROLS

- 1500M6 Prop RIM Head Control w/ MK6i Display (E2327) . 2,875.00
- SMHD Small HeadrestNo Charge
- LGHD Large HeadrestNo Charge
- 1558MM6 Compact Joystick w/Display (E2373) 1,130.00

ALTERNATIVE DIGITAL CONTROLS

- SNPM6 Sip-N-Puff Digital Interface Box (E2325) 1,500.00 ²¹
- 5018M6 Wafer Board & Digital Interface Box (E2322) 1,550.00 ²¹
- 5020M6 Mini Tash Joystick & Digital
Interface Box (E2321) 1,750.00 ²¹

SIP-N-PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit (E2326) .380.00

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick – Less Std Mount 770.00
- 5018 Wafer Board..... 310.00 ²²
- 5020 Mini Tash Joystick..... 405.00 ²²
- PACM6 Prop. Attendant Control w/ quick release mount. 850.00 ²²
- 1552M Digital Attendant Control 850.00

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt MPJ+ (E1028)..... 279.00 ²³
- QLA Quad Link Retract. Joystick Mnt SPJ+ (E1028)..... 279.00 ²³
- ARM250 Stealth Height Adj. Joystick Mnt (E1028) 210.00
- LEFT Left Handed Mounted JoystickNo Charge
- RIMHW RIM Hardware For 1500 (E1028) 219.00
- PKG32669 Complete Bib Assembly (Tash/1558) 265.00
- GATMPJM6 Gatlin Midline for MPJ (E1028) 600.00
- GAT1812 Gatlin Midline for 1812 (E1028)..... 600.00
- RCM Rear Cane Mount for Compact Joystick 395.00

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer w/Pro Memory Card 475.00
- 1813SD Pro Memory Card (USB Ready) 90.00
- 1560 "T" Handle Flexible Joystick Ext. (E2323) 77.00
- 1561 Straight Handle Flexible Joystick Ext. (E2323)..... 77.00
- 1826 Chin Cup (2324) 48.00
- A24VPS 24V MK6 Auxiliary Power Source 200.00
- AUX12M6 Auxiliary Module For 1 & 2..... 795.00
- AUX34M6 Auxiliary Module For 3 & 4..... 795.00
- MPY Monoport "Y" Cable for MPJ+ or Display 15.00
- 2PB Dual Push Button..... 147.00 ²⁴
- 2SW Dual Toggle 172.00 ²³
- PBMT Push Button Mounting Kit 252.00
- IR Infrared Module with Wireless Mouse Control 995.00 ²⁵

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch – Black 156.00
- ASL304 Light Wobble Switch w/611 Mount 432.00 ²⁶
- ASL202J Fiber Optic Switch..... 610.00 ²⁶
- ASL208J Adjustable Proximity Switch..... 550.00 ²⁵
- ASL504M6 Remote Emergency Stop Switch 1,185.00
- ASL611 Adjustable Switch Mounting Bracket 252.00

FOOTNOTES

- ²¹ Must also use either a Multiple Drive Joystick or MK6i Display.
- ²² Requires SNPM6 Sip-N-Puff Interface Box.
- ²³ QLAM6 and QLA come standard with ARM250. Not available with TRRO.
- ²⁴ Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options
- ²⁵ Capable of controlling up to six remote Infra-Red controlled devices and also Computer Mouse operation through the driver control. Includes stereo switch (with hook/loop mount) for mouse clicks
- ²⁶ Not available with SPJ+ Joystick. Must add MK6 Auxiliary Power Source.

NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.

Invacare® TDX® SI Rehab Seat Power Wheelchair Accessories and Options

MISCELLANEOUS

BACK ACCESSORIES

- L8601 Omit Back Upholstery No Charge
- RBP Reinforced Back Upholstery No Charge

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks 45.00
- SEAT02 Oxygen Holder (ADJASBA Seat Only) (E2208)..... 150.00
- ADD Adductors Small (E0956)..... 200.00
 - ADDHW Hardware for ADD (E1028)..... 400.00
- ADD L Adductors Large (E0956)..... 200.00
 - ADDHWL Hardware for ADD L (E1028)..... 400.00

BATTERY TRAY TYPE ²⁷

- 22TRY 22NF Style Tray **Standard**

INSTALLATION OF BATTERIES

- 22NFBATTERY 22NF Battery (On Chair) (E2361).....300.00

BATTERY CHARGER OPTION

- 110CHARGER 110 Volt Battery Charger No Charge
- LESTR Lester 8 Amp Battery Charger (E2367) 320.00
- L8900 Omit Battery Charger Less 80.00

VENT TRAY

- COMP Compact Vent Tray (E1030)..... 1,540.00 ²⁸

FRONT RIGGINGS

HANGER BRACKET TYPE

- ASB Center Pivot Style Bracket No Charge
- PSB Pin Style Bracket No Charge

FRONT RIGGING ²⁹

CENTER MOUNT

- CTMT-F Center Mount Adj. Knee Angle 300.00 ³⁰
- CTMT-FJ Junior Center Mount Adj. Knee Angle 300.00 ³¹
- SM Small Footplates (5.25" x 7.5").....No Charge
- MED Medium Footplates (5.50" x 9.5").....No Charge
- LRG Large Footplates (5.90" x 11.5").....No Charge
- CALF Calf Pads For Center Mount.....200.00

CENTER PIVOT STYLE (PAIR)

- 70NHD 70° Swingaway Footrests 13.5"-17.5" **Standard** ³²
- 70STAPER 70° Tapered Footrests 13"-17"..... No Charge ³³
- AT5544 Manual ELR (E0990)..... 255.00 ³⁴

PIN STYLE HEMI (PAIR)

- PHW93 70° Swingaway Footrests 13.5" - 19".....No Charge ³⁵
- PHAL4A Manual Articulating/ELR (K0053)..... 225.00 ³⁶

FOOTNOTES

²⁷ 22TRY Base Width 24".

²⁸ Approximate dimensions are 11" W x 4" D x 13" H. Designed for use with the latest and most compact style vent systems.

²⁹ All front rigging options include a standard footplate.

³⁰ Allows 60 to 90 degree static knee angle positioning with a drop back lever for easier transferring. Extension range 12.5" – 15". Must choose either SM, MED or LRG footplate.

³¹ Allows 60 to 90 degree static knee angle positioning with a drop back lever for easier transferring. Extension range 9.25" – 12". Must choose either SM, MED or LRG footplate.

³² Standard with the 1651 composite footplate. With seat widths of 21" to 22" and depth of 16", must pick footplate 1350. The AT5543 is not available.

³³ Standard with impact guards and adjustable angle flip up footplates (AT5543).

³⁴ Not available with 16" wide by 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.

³⁵ Standard with 1651 composite footplate.

³⁶ Not available on 12" or 13" width seats.

FOOTBOARD ³⁷

- 90FB Flip up removable footboard 150.00

FOOTPLATE OPTIONS (PAIR)

- 1651 Composite Flip Up Footplate No Charge ³⁸
- 1350 Extra Large Aluminum Footplates No Charge ³⁸
- AT5543 Adjustable Angle Flip-Up Footplates (K0040).... 165.00 ³⁹
- CTMTFB Footboard for Center Mount..... 150.00 ³⁹

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap (K0038) 27.00
- 1600BK Heel Loops w/Ankle Straps (pair) (E0951)..... 56.00
- IMPACTG Impact Guards..... 51.00
- LPT2 Longer Pivot Slide Tube for 1350..... 83.00
- ALPT Longer Pivot And Slide Tube..... 84.00
- ALPT4 Longer Pivot And Slide Tube 4" 84.00

FOOTNOTES

³⁷ Footboard only available with SEMIRB seat style.

³⁸ Not available with AT5544 and 16" wide by 16" – 18" deep or 17" wide by 16" deep and low seat-to-floor.

³⁹ One piece footboard for Center Mount (CTMT-F) Riggings.

NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.

Invacare® TDX® SI Rehab Seat Power Wheelchair Accessories and Options

TDX SI Specifications		
Seat Width Range: (In 1" increments)		
	Adult (Adjust)	Junior (Adjust)
ADJASBA	16" – 24"	12" – 16"
SEMIRB (VAN)	16", 18", 20", 22"	
Seat Depth Range:		
ADJASBA	12" – 22"	12" – 16"
SEMIRB (VAN)	16" – 20"	
Back Height Range:		
	12" – 24"	12"-22" Adj. Height
	18", 20", 22", 24", 26"	Recliner Back Heights
	25" & 29.5"	Contoura Back
	18" w/o headrest	SEMIRB
	25" w/headrest	SEMIRB
Back Angle Range:		
	85° - 105° in 5° increments * except Recliner	
Seat-to-Floor (approximate):		
<u>Standard:</u>	16.5" @ 0°	
<u>Medium:</u>	18.5" @ 0°	
<u>Tall:</u>	20.5" @ 0°	
	At 5° add .5"	
Overall Width:		
	24" (w/o Joystick)	
Overall Height (with 20" back height & 0° seat angle)		
<u>Standard:</u>	35.5"	
<u>Medium:</u>	37.5"	
<u>Tall:</u>	39.5"	
	38" w/o headrest for SEMIRB	
	45" w/headrest for SEMIRB	
Overall Length:		
	35.25" w/o front rigging	
Weight:		
<u>2 Pole Motor</u>		
<u>ADJASBA</u>	w/o 22NF Batteries:	166 lb
	With two (2) 22NF Batteries:	264 lb
<u>SEMIRB</u>	w/o 22NF Batteries:	179 lb
	With two (2) 22NF Batteries:	277 lb
<u>Drive Wheels/Tires:</u>	(Foam Filled Standard)	14"x 3"

MEMORY CARD – BASIC

- Standard with ALL expandable control systems (MPJ+, PSR+, PSF+, MK6i Display).
- Used to back up/restore programmed settings for only one chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not contain advanced diagnostics/help screens/file structure
- Not compatible with SPJ+ Joysticks.

MEMORY CARD – PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, help screens, software updates, file storage/retrieval.
- Used to back up multiple chairs and programming settings.
- Intended for used by qualified providers only.
- Not compatible with SPJ+ Joysticks.

Batteries
(not included)

22NF batteries with terminal configuration (positive on the left and negative on the right) MUST be used. 22NF batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection.

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
* Not available with Contoura Seat Cushion							
Back Heights 25" or 28"							

Weight Limit 300

Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the products.

NOTE:
All weight limits are inclusive of add-on equipment and seating systems. All specifications and dimensions are approximate.

SEATING AND POSITIONING

**FOR YOUR SEATING NEEDS
PLEASE REFER TO THE FOLLOWING:
REHAB SEATING AND POSITIONING PRODUCTS
SEE PRICE LIST FORM 07-027.**

Cantilever Arm Height Adjustment			
Back Type	Height	Adjustment	Adjustment with laterals
CURV	20" or greater	11" to 16"	11" to 14"
CURV	19"	11" to 15"	11" to 13"
CURV	18"	11" to 14"	11" to 12"
CURV	17"	11" to 13"	no adjustment
CURV	16"	11" to 12"	not available
Contoura	25" & 28"	12.5" to 14.5"	no adjustment

Invacare Corporation www.invacare.com



USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Products designated by ASL are Manufactured by Adaptive Switch Labs, Inc. Therafin is a registered trademark of W/R Medical Electronics Co. Otto Bock is a trademark of Otto Bock Orthopaedische Industries GMBH & Co. Specifications and prices are subject to change without notification.

©2008 Invacare Corporation